

Product Name : Chocolate Fountain Machine	Product Code : NLAB-VOCATIONALTVET270103
-----------------------------------------------------	----------------------------------------------------

--

Description :

Chocolate Fountain Machine

Technical Specification :

Capacity : 2 lb
Material : Stainless Steel
Voltage : 210 V
Frequency : 60 Hz
Power : 2.5 kW

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

