

Product Name :
Cholera Bed Foldable

Product Code :
NLBMEDICAL-M2660001



Description :

Cholera Bed Foldable

Technical Specification :

Cholera Bed Frame is made of Mild Steel Tube.

The base of the Cholera bed is made of Nylon Backing With Provision of Hole in the Middle.

It is foldable and light in weight .

There is a provision of IV rod.

It comes with Plastic Sanitary Bucket.

The carrying capacity of the bed is 120 kg.

Product Size : 180 L x 60 W x 55 H cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Cholera Bed Foldable at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003