

Product Name :
Chromatography Column with sintered disc and P.T.F.E. key stopcock

Product Code :
NLAB-SCIENTIFICEQ18006



Description :

Chromatography Column with sintered disc and P.T.F.E. key stopcock

Technical Specification :

Plain with P.T.F.E. screw type Rotaflow stopcock

Effective Length

- ~~100~~ Non
- ~~300~~ 300.01
- ~~200~~ 200.02
- ~~300~~ 300.03
- ~~200~~ 200.04
- ~~300~~ 300.05
- ~~500~~ 500.06
- ~~800~~ 800.07
- ~~800~~ 800.08

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments