

Product Name : Chromatography Column with sintered disc and Screw Type P.T.F.E. Rotaflow stopcock	Product Code : NLAB-SCIENTIFICEQ18004
---	---




Description :

Chromatography Column with sintered disc and Screw Type P.T.F.E. Rotaflow stopcock

Technical Specification :

Chromatography Column with sintered disc and Screw Type P.T.F.E. Rotaflow stopcock.

Effluent Length mm
300.01
300.02
300.03
300.04
300.05
500.06
800.07
800.08



Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments