

Product Name :
Circular Saw and Accessories

Product Code :
NLAB-MACHINELABWK34011



Description :

Circular Saw and Accessories

Technical Specification :

Type: Cold Saw for metal cutting

Cutting Saw (type): Circular saw

(Size): Diameter – 304.8mm

Capacity (Steel/90 degrees): 51mm

Speed: Approx. 44 rpm

Head Rotation: +/-45o

Pedestal Stand: Must be well equipped

Features: Mechanical coolant pump

Adjustable blade stop

Power Source: 220V 50Hz, Single Phase

Accessories:

Standard Accessories – including Tool Kit 3 set

Spare Circular Saw – Dia. 304.8 (5 pcs)

Operations manual

Standard Welding Equipment Laboratory Set Circular Saw and Accessories
Manufacturer, Supplier and Circular Saw and Accessories Exporter in India

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003