

Product Name :
Circulating Water Bath with Cooler Unit

Product Code :
NLAB-CEMETLAB130037



Description :

Circulating Water Bath with Cooler Unit

Technical Specification :

Automatically maintains the set temperature

Two channels electronic thermoregulator: heating/cooling with digital display

Complete with recirculating unit

Incorporating refrigeration compressor

Ideal for conditioning Asphalt specimens for Marshall and Indirect tensile tests

Ideal for cement specimens cure in water

A multi-purpose unit usable for many other applications

Capacity: 40 l

Power: 2000 W

Temperature range: +5 to +60°C

Accuracy: $\pm 1\%$

Inside dimensions: 550x360x200 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003