

Product Name :
Clamping Machine

Product Code :
NLAB-MACHINELABWK33026



Description :

Clamping Machine

Technical Specification :

(30) Total 3-1/2" High jaw 40" Opening T-Handle Clampss

(5) Total Gluing Levels

Free standing Unit with clamps hanging at a 20 degree down angle

Vertical Clamp Hanger

100% Steel Construction for long lasting durability

Additional Clamp Arms available upon request - Maximum of (50) clamps

Wood Workshop Portable Tools Clamping Machine Manufacturer,Supplier and Clamping Machine Exporter in India

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003