

Product Name :
Clinical Lab Rotary Shakers

Product Code :
NLABMEDICAL-M3600001



Description :

Rotary Shakers

Technical Specification :

Clinical Lab Rotary Shakers - These Shakers are ideal for mixing and development of cultures, chemicals, solvents, and assays in Microbiological, Cell Culture and Life Science laboratories, The Unit has base assembly fabricated from heavy mild steel sections. Brushless Induction drive motor with frequency drive makes the unit suitable for continuous non-stop operation. Step less electronic frequency control ensures gentle start and maintains preset speed. Compact counter balanced drive mechanism ensures high stability and reliability even in continuous operation & uneven load distribution.

4? Large LCD display

Universal Platform with rounded corners to accommodate interchangeable clamps of assorted sizes

Brushless induction motor with VFD speed control

Elegant Plastic front

Variable speed 20 -300 RPM

Shaking amplitude 40 mm

Auto start after power resumption

Platform Size 24" X 24"

Max. Volume 18 Litres

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003