

Product Name :
Clip On Meter Clamp Meter

Product Code :
NLAB-TESTINGLABM18002



Description :

Clip On Meter Clamp Meter

Technical Specification :

Technical Specifications :

Measurement : ACA, DCA, ACV, DCV, ohms, Diode

DC voltage range : 200mV ~ 600V

Resolution : 0.1mV ~ 1 V

AC Voltage Range : 200V ~ 600V

Resolution : 0.1mV ~ 1V

Resistance : 2000 ohms

Resolution : 0.1A ~ 1A

Built in data hold function

Meet IEC 1010 safety measurement

Built in low battery indicator

High internal impedance

Overload protection circuit for all range

LSI circuit for high reliability and durability

Current sensor : Hall effect sensor

Zero adjustment : Automatic

Sampling time : 0.4 sec.

Power : DC 9 Volt battery

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003