


Product Name :
Clipping and Clamping Using PN Diode Trainer

Product Code :
NLAB-ELECTRICALAB210028



Description :

Clipping and Clamping Using PN Diode Trainer

Technical Specification :

SCOPE OF LEARNING:

POSITIVE CLIPPER

NEGATIVE CLIPPER

POSITIVE CLAMPER

NEGATIVE CLAMPER

TECHNICAL SPECIFICATIONS:

Power Supplies:

DC Supply 0-5V DC, 150mA.

Operated on Mains power 230V, 50Hz +10%

Components are mounted on the panels are:

RESISTANCE

Diodes

Capacitors

Voltage Control through Potentiometer.

Series Resistor 470E

SALIENT FEATURES:

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Instruction manual.

Connections are brought out through 4mm Colored Sockets.

Patch Cords 4mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003