

**Product Name :**  
CNC Milling Machine

**Product Code :**  
NLAB-MECHANICALAB110002



**Description :**

CNC Milling Machine

**Technical Specification :**

NaugraLabEquipments is a leading educational cnc milling machine manufacturers, suppliers and exporters company in India. We have been supplying an extensive variety of educational cnc milling machine around all over the world at the very affordable cost.

Our CNC lab milling machine are designed and manufactured by a dedicated team of professionals to provide superior, high-performing and long lasting service.

**Specification**

units:x6140 CNC Router, 3-axis setup, Windows 7 or latest operating platform, 19" flat screen monitor with 2 USB connections. Easy input of processing data via elucam 3-D display

Table size: 1600mm x 400mm

Number of T slots:3

Maximum load:800kg

Longitudinal travel:(both manual and auto)900/800

Cutting feed speed: mm/min 19-950

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003