

Product Name :
CNC Milling Trainer Deluxe

Product Code :
NLAB-MECHANICALAB150005



Description :

CNC Milling Trainer Deluxe

Technical Specification :

The Deluxe LabTek-PT-10005 metric comes equipped with all standard features upgraded to include a laser engraved 12" (305mm) solid aluminum base, laser engraved scales on the "X" axis table, 2" (51mm) hand wheels on the "X" and "Y" axes, a 2- 1/2" (63mm) adjustable "zero" handwheel with ball bearing thrust on the "Z" axis, a mill head stock spacer to provide more throat distance, and drill chuck and drawbar.

Specification:

8100" (2060mm) table to spindle :

2125" (5400mm) (without headstock spacer) :

Fitted with 50# (9070g) spacer block:)

865# (220# Axis) :

500# (127# Axis) :

625# (159# Axis) :

405" through spindle :

Spindle thread :

Spindle taper :

Blade wheel graduations :

Width (476mm)

Depth (476mm)

Height (560mm) :

Table size (70x330mm)

Motor provision :

36 lb. (16.3 kg) IMovements in addition to

Shipping weight :

Headstock rotation 90° left/right

Movements in addition to X-, Y- and Z-axes :

90 Volt DC with electronic speed control

Motor :

220VAC, 50/60 Hz

input:

70-2800 RPM continuously variable by electronic speed control

Spindle speed range :

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003