

**Product Name :**  
Coarse Aggregate Density Test Set

**Product Code :**  
NLAB-VOCATIONALTVET17007



**Description :**

Coarse Aggregate Density Test Set

**Technical Specification :**

Buoyancy Balance, of capacity 8000 g for aggregates between

63mm and 5mm

8" Brass Sieve no.4 Mesh

1 Pan, Galvanized

Power 220 vAC, 50/60H

Density Basket: with dimensions 203 mm diameter x 203 mm deep.

constructed of reinforced no. 8 wire mesh.

Suspension hook

Water tank with dimensions: 400 x 400 x 330 mm approx. (L x W x H)

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003