

Product Name :
Coated Fabric Testing Lab Equipment

Product Code :
NLAB-TESTINGLABM50001



Description :

Coated Fabric Testing Lab Equipment

Technical Specification :

Size of the Test Specimen: 40 X 15 mm

Creasing Load: 1.0 kg

Area of Loading: 15 X 15 mm

Distance between edge of the

Clamp to the axis: 2 mm

Crease Recovery Angle: 10 °C â€“ 180 °C X 0.5

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

