

Product Name :
Coaxial Binocular Microscope

Product Code :
NLAB-EDUSCHOOLEQ173072

Description :

Coaxial Binocular Microscope

Technical Specification :

Rectangular 115 x 125 mm with co-axial mechanical stage.

Box type base with 6v-20w, variable halogen illumination.

Sub stage movable by rack and pinion Abbe Condenser of 1.25 N.A. with iris diaphragm and filter

Achromatic objectives: 4x, 10x, 40x, 100x(SL).

High quality Binocular Microscope for School, Colleges and Research Institutes.

Eyepieces: 10x(pair) WF.

Huygenian eyepieces: 5x (pair).

Provided with a Binocular observation head, inclined at 45 degree, with interpupillary distance 55mm to 75 mm & diopter adjustments.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003