

|  |   |
|--|---|
| <b>Product Name :</b><br>Coaxial Binocular Microscope  | <b>Product Code :</b><br>NLAB-EDUSCHOOLEQ173072 |
|  |   |
| <b>Description :</b>   |   |
| Coaxial Binocular Microscope   |   |
| <b>Technical Specification :</b>   |   |
| Rectangular 115 x 125 mm with co-axial mechanical stage.   |   |
| Box type base with 6v-20w, variable halogen illumination.  |   |
| Sub stage movable by rack and pinion Abbe Condenser of 1.25 N.A. with iris diaphragm and filter                                      |   |
| Achromatic objectives: 4x, 10x, 40x, 100x(SL).   |   |
| High quality Binocular Microscope for School, Colleges and Research Institutes.  |   |
| Eyepieces: 10x(pair) WF.   |   |
| Huygenian eyepieces: 5x (pair).  |   |
| Provided with a Binocular observation head, inclined at 45 degree, with interpupillary distance 55mm to 75 mm & diopter adjustments. |   |

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003