

Product Name :
Coaxial Trinocular Microscope

Product Code :
NLAB-EDUSCHOOLEQ1730147

Description :

Coaxial Trinocular Microscope

Technical Specification :

Illumination: Fitted with built in base halogen light source, fitted with 6 v-20 W Halogen bulbs, On / Off Switch and light intensity regulator, complete with power cord and plug. Also a set of plano-concave reflector provided for day light use.

Sub Stage: Abbe N.A.1.25 condenser, focusable with rack and pinion through 25 mm continuously visible iris diaphragm with built in swing out filter holder.

Head: Inclined Trinocular head

Body: Inclined Binocular Body, Rotatable through 360o

Focusing: Single Co-axial knob provided based on a four step gear reduction system for coarse and fine focusing for silk smooth operation.

Mechanical Stage : Large stable base (120x 125mm) fitted with co-axial precision graduated mechanical stage for X and Y movement of slide upto 50mm x 75mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003