

Product Name :
Cold Box 22.40 Ltr

Product Code :
NLABMEDICAL-M2330001



Description :

Cold Box 22.40 Ltr

Technical Specification :

Cold Box 22.40 Ltr

Cold Box 22.40 Ltr external surface material is made of LLDPE

Cold Box 22.40 Ltr internal lining material is made of LLDPE

Cold Box 22.40 Ltr is also known as Large Cold Box, Long Range

WHO Reference PQS Number: E004/046

The Vaccine Storage Capacity is 22.40 Litres

The fully loaded weight is 49 Kgs.

The Empty Weight with empty icepacks is 21.55 Kgs.

The Insulation Material is made of CFC free Polyurethane

The thickness of Insulation is 100 mm

The External Dimensions are (L x W x H) 77 x 61.5 x 51.5 cms

The Internal dimensions are (L x W x H) 52 x 37 x 28.2 cms

The Vaccine Storage Dimensions are 44.8 x 29.8 x 16.8 cms

The storage box and the Lids are Hinged & Removable

The total Number of Ice packs required are 31 x 0.6 Litre which is supplied with the unit

Ice-packs type: WHO PQS No. E005/009

The Cold Life Requirement as per WHO at 43 deg C \pm 1 deg C without opening should be Minimum 153.5 Hours which satisfies this model.

Warm Life at -20oC is 69 Hours

Cool Life at +43oC is 41.7 Hours

We are Manufacturer, Exporters, and Suppliers of Cold Box 22.40 ltr at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003