

**Product Name :**  
Cold Box Long Range

**Product Code :**  
NLABMEDICAL-M2300001



**Description :**

Cold Box Long Range

**Technical Specification :**

Type: Cold Box, Long Range

Vaccine Storage Capacity As per Standard PQS load: 6.0 Litres

Gross Storage Capacity: 11.61 Litres

Weight-Fully Loaded: 25,40kg

Weight-Empty (with empty Ice Packs): 12,82kg

External Surface Material: LLDPE (Linear Low Density Polyethylene)

Internal Lining Material: LLDPE (Linear Low Density Polyethylene)

Insulation Material: CFC Free Polyurethane

Insulation Thickness: 100mm

External Dimension (L X W XH): 62,40 x 50,20 x 42,60cm

Internal Dimension (L X W XH): 40,30 x 28,20 x 21,10cm

Vaccine Storage Dimensions: 33,00 x 21,20 x 16,60cm

Lid Type & Fixings: Hinged

Number of Ice-packs Required: 26

Number of Ice-packs Supplied: 26

Ice-pack type: 0.4L compliant with PQS E0S/IP0I.1

Volume of Ice-pack: 0.4 Litres

Performance: Min. Cold Life Std. Req. 30hrs

Cold life without opening at +43C: 132 hrs 18 minutes

Warm Life at -20 C: 48 hrs 41 minutes

Cool Life at +43 C: 29 hrs 47 minutes

Shipping Dimensions

Gross Weight (Approx.): 16kg

Size of carton (LxWxH) (approx.): 65 x 53 x 46cm

Volume of each carton (approx.): 0.16 CBM

No. of Units/20 Container: 175

No. of Units/40 Container: 375

No. of Units/40 High cube Container: 445

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003