

<b>Product Name :</b> Cold Room	<b>Product Code :</b> NLABMEDICAL-M3180005
<b>Description :</b>	
Cold Room India, Cold Room manufacturer India, suppliers, Medical Cold Room India, Hospital Cold Room India	
<b>Technical Specification :</b>	
Cold room has total body with heavy sturdy angle iron & channels outer body is made of MS sheet duly painted.  Inner chamber made S.S. single door opening system with door opening alarm.  Temperature 4° C is brought down by hermetically sealed compressor and controlled by microprocessor based temperature variation alarm facility controller, puff insulation to maintain the temperature, inner chamber illumination total unit is moveable on casters.	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003