

Product Name :
Color Fastness Tester Machines

Product Code :
NLAB-TESTINGLABM59001



Description :

Color Fastness Tester Machines

Technical Specification :

Features

- Imported Mercury-Tungsten Lamp used
- Economical General Duration testing machine for sunlight
- Mercury tungsten light for simulation of light fading by sunlight
- 12 sample holders included can be tested at same time and 10 samples can be packed in each holder.
- Equipped with Rapid fade unit.
- Exposure Area : 1800cm²
- Lamp Life: 5000h

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003