

Product Name : Color Scanner Trainer	Product Code : NLAB-VOCATIONALTVET31006
--	---

--

Description :

Color Scanner Trainer

Technical Specification :

Each trainer is designed to instruct service technicians in the operation, theory, servicing, problem diagnosis and repair of networking and communications devices for today computer I.T. environment.

All of the trainers are based upon hardware machines which have been modified to facilitate training and have a special electronic fault insertion system.

This trainer is part of I.T. training system.

Naugralabequipments
Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003