

Product Name :
Color Television Trainer, 21" Color

Product Code :
NLAB-EDUSCHOOLEQ194002



Description :

Color Television Trainer, 21" Color

Technical Specification :

These trainers can be used as stand-alone units or combined with our VCR trainers to make a complete VCR/TV training system.

Special attention is given to fault diagnosis to facilitate the teaching of trouble shooting techniques.

They are specifically created to teach television operation, theory, repair and servicing

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003