

Product Name :
Colostomy Bag

Product Code :
NLABMEDICAL-M3040001



Description :

Colostomy Bag

Technical Specification :

- Colostomy Bag is used for closed colostomy procedures.
- Colostomy Bag is manufactured from clinical grade PVC sheet.
- Colostomy Bag is coated with Hypo-allergic adhesive to prevent skin maceration.
- Colostomy pouching system is a medical device prosthetic that provides a means for the collection of waste from a surgically diverted biological system (colon, ileum, urinary) and the creation of a stoma.
- Pouching systems are most commonly associated with colostomies, ileostomies, and urostomies.
- Colostomy Bag is Non Sterile.

Packing

- Pack of 10 No
- Master Carton of 1000 No

- Master Gross Weight: 13 Kg
- Master Carton Size: 58 x 16 x 34 cm

We are Exporters and Suppliers of Colostomy Bag at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003