


**Product Name :**  
Commercial Refrigeration Trainer

**Product Code :**  
NLAB-TECHNICALAB11004



**Description :**

Commercial Refrigeration Trainer

**Technical Specification :**

- Overall size: 100"W x 27"D x 74"H
- Electrical requirements: Single phase, 220V AC, 60 Hz, 40A (115 VAC, 60 Hz for water tower)
- Utilities required are city water, drain to water, and means of venting water tower
- Plumbing to city water and water tower not included
- Weight: 1300 lbs.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

