

Product Name :
Compact Infrared Thermometer Digital

Product Code :
NLAB-EDUSCHOOLEQ208001



Description :

Compact Infrared Thermometer Digital

Technical Specification :

IR-thermometer with laser enables fast, easy and accurate readings for surface temperature measurements. With built in laser pointer for non-contact infrared temperature measurements of hot or moving objects from a safe distance or difficult measuring points. With this infrared thermometer, it is also possible to find leaks or cold bridging, for example in a wall or. The automatic color change function allows quick and intuitive assessment and indicates where must be sealed.

Measuring range -50...+380 °C

Accuracy +/- 2.0% + 2 °C

3 1/2 digit, 15 mm multifunction LCD-display with backlight and data hold

Unique flat surface, modern housing design

Max. indication: 1999, overload-indication and beeper

Fast temperature difference measurement with LED indicator (red, blue or green LED)

°C/°F switchable

Built-in laser pointer identifies target area (class II; <1 mW)

Overrange indication

Auto data hold, auto power off

Safety: IEC-60825-1

Accessories: carrying case, battery and manual

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003