

Product Name :
Compressible Flow Bench

Product Code :
NLAB-VOCATIONALTVET19002



Description :

Compressible Flow Bench

Technical Specification :

Compressible Flow Bench

Pressure measurement bench allows the student to investigate the measurement of one of the fundamental parameters that will be present in almost every branch of engineering and physics.

The unit should be bench top and self-contained having its own means of pressure generation.

Panel mounted U tube and inclined tube manometer together with a positive pressure bourdon tube pressure gauge and compound pressure

gauge.

The manometers should allow investigation of the use of U and inclined tubes for pressure measurement and demonstrate the use of fluids of different density.

All of the manometers and the panel mounted pressure gauges may be interconnected and linked to the common pressure source supplied

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003