


Product Name :
Compression Balance (Plastic)

Product Code :
NLAB-SCIENTIFICEQ11006



Description :

Compression Balance (Plastic)

Technical Specification :

- Economical top loading weighing scales having small footprint, and meant for quick, general purpose weighing with sufficient accuracy.
- Slanted wide view circular dial and large pointer needle provides error free and convenient reading of scales even from a distance.
- Molded plastic body in beautiful pastel colors and trendy matching dials and large scoop pans at the top makes its use convenient for weighing small parts, powdered or granulated items.
- With zero adjustment knob.
- **Capacities available are:**
- 500gx2g

- 1kgx5g
- 2kgx10g
- 5kgx25g

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003