

Product Name :
Compression-Flexure Cement Testing Frames

Product Code :
NLAB-CEMETLAB130039



Description :

Compression-Flexure Cement Testing Frames

Technical Specification :

Single compression station, 500 kN or 600 kN cap.: for compression tests on low strength specimens.

Double station model, 600/15 kN or 500/15 kN cap.: extends the possible applications to flexural tests on cement and compression tests on low strength specimens.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003