

Product Name : Compression Testing Equipment Electrically cum Manually (Hand) Operated	Product Code : NLAB -CONCRLAB120007
--	---



Description :

Compression Testing Equipment Electrically cum Manually (Hand) Operated

Technical Specification :

Compression Testing Equipment Electrically cum Manually (Hand) Operated

Capacity (In kN)			
Pressure Gauge (In kN)			
Piston Count (In kN)			

Pressure Gauge Diameter (mm)			
Stroke (mm) Daylight (mm)			
Horizontal Clearance (mm)			
Stroke Diameter (mm)			
Stroke Diameter (mm)			
Stroke Travel (mm)			
Motor H.P			
Motor Phase 50Hz			
Specimen Size (Can be Tested):			

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003