

**Product Name :**  
Compression Testing Equipment Manually Hand Operated

**Product Code :**  
NLAB -CONCRLAB120006



**Description :**

Compression Testing Equipment Manually Hand Operated

**Technical Specification :**

Capacity (In kN): 1000 kN

Pressure Gauge (In kN): 1000 kN

Least Count (In kN): 5 kN

Pressure Gauge Diameter (mm): 150

Vertical Daylight (mm): 320 (Adj)

Horizontal Clearance (mm): 260

Platen Size (mm): 220 x 165

Ram Diameter (mm): 165

Ram Travel (mm): 50

Cube: 100mm, 150 mm

Cylindrical: 100 mm x 200 mm, 150 mm x 300 mm

Brick: 100 x 100 x 220 mm

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003