

Product Name :
Compression Testing Machine 2 Pillars Type

Product Code :
NLAB-COMPSLAB280006



Description :

Compression Testing Machine 2 Pillars Type

Technical Specification :

Compression Testing Machine – 2Pillars Nut Type Model, complete with Single Pressure Gauge & Self Made Motor, 440 'V' Load Frames: Available in 100 Tons, 150 Tons, 200 Tons. Model Available: Hand Operated & Electrically cum Hand Operated

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003