

Product Name :
Compression Testing Machine Channel Type

Product Code :
NLAB-COMPSLAB280008



Description :

Compression Testing Machine Channel Type

Technical Specification :

Compression testing machine – channel type model, complete with single pressure gauge & self made motor, 440 'v' load frames: available in 100 tons, 150 tons, 200 tons. Models available: hand operated & electric cum hand operated.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003