

Product Name :
 Compression Testing Machine for Concrete

Product Code :
 NLAB-ENGINEERINGLB33002



Description :

Compression Testing Machine for Concrete

Technical Specification :

Compression Testing Machine for Concrete Lab Equipments Manufacturer, Suppliers & Exporters Compression Testing Machine For Concrete .

(2 Pillar Model): Compression Testing Equipment Manually (Hand) Operated fitted with single load gauge. The loads are measured on Bourdon tube type load gauges which are calibrated against certified proving rings. The load gauges are fitted with a maximum load pointer

Capacity (in KN)

1000 KN

1000 KN Gauge (In KN)

5000	Load Count (In KN)
200	Pressure Gauge
	Diameter (mm)
3/16"	Vertical Daylight (mm)
250	Horizontal Clearance(mm)
240	Core Size (mm)
165	Pen Diameter (mm)
50	Pen Travel (mm)
Specimen Size (Can be Tested)	
100	100mm x 150mm
100	Cylindrical 200mm
	150mm x 300mm
100	Block 100 x225 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003