

**Product Name :**  
Compressor Analyzer

**Product Code :**  
NLAB-VOCATIONALTVET13013



**Description :**

Compressor Analyzer

**Technical Specification :**

Tied-up compressors may be freed by reversing motor rotations

Tests windings for opens, continuity, and grounds

Tests start and run capacitors and system relays by substitution

Measures running amps, 0-10, 50A range

Measures start and run volts, 0-150, 600V AC ranges

Compressors are tested in the manner in which they were designed to operate: Split-phase, capacitor start, or PSC

Has provisions for connecting external run capacitor.

30A breaker protected to allow high starting current with adequate overload protection

Complete with 220V adapter cord and capacitor leads

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003