

Product Name :
Compressor Controls Trainer

Product Code :
NLAB-ENGINEERINGLB46001



Description :

Compressor Controls Trainer

Technical Specification :

Compressor Controls Trainer Lab Equipments Manufacturer, Suppliers & Exporters Technical Description
The benchtop unit was designed for electrical training for refrigeration system installers and is particularly suitable for exercises in smaller groups. The electrical connection of a commercially available refrigerant compressor using a circuit diagram is covered. The starting, protection and control devices necessary for safe operation are arranged on the front panel of the unit. Test points provide convenient access for measuring the electrical circuit.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003