

<b>Product Name :</b> Computer Compatible Manometer Bank	<b>Product Code :</b> NLAB-TECHNICALAB26005
<b>Description :</b> Computer Compatible Manometer Bank	
<b>Technical Specification :</b> Electronic multiple channel pressure measurement system designed to replace liquid manometer banks. 16 simultaneous measurements displayed on a PC (not supplied). Sophisticated software included, for logging, recording and displaying data. 16 channels of 350mm H2O capacity, or 14 channels of 350mm plus 2 channels of 125mm H2O.	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003