

Product Name :
Computer Embroidery Machine

Product Code :
NLAB-EDUCATIONTVET21012



Description :

Computer Embroidery Machine

Technical Specification :

Computer: 8" LCD touch screen panel

Main motor: Dahao servo motor

Auto trimming, Automatic color changing and start up

Automatic thread break detection

Stitch length: 0.1-12.7mm

Storage capacity: single head - 16,000,000 stitches and 400 designs

Multi-language operating interface

Working Speed: 400~1200rpm

Rotary hook: Japanese Hook

Power: 110V / 220V, 50HZ / 60HZ

Upper thread holding device

head single

Embroidery area : 600x400 mm

Embroidery shirt area : 300x300mm

Power 1600 watt

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003