

Product Name :
Computer Networking and Communication Trainer

Product Code :
APPNLAB-
ENGINEERINGLB54054



Description :

Computer Networking and Communication Trainer

Technical Specification :

Each trainer is designed to instruct service technicians in the operation, theory, servicing, problem diagnosis and repair of networking and communications devices for today computer I.T. environment.

All of the trainers are based upon hardware machines which have been modified to facilitate training and have a special electronic fault insertion system.

This trainer is part of I.T. training system.

We are well-known manufacturers, OEM suppliers of Computer Networking and Communication Trainer for Applied Machines. Contact us for high quality Computer Networking and Communication Trainer for Applied Machines for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003