

Product Name :
Concrete Cube Mould

Product Code :
NUG -CTLE0011NLAB
-CONCRLAB120011



Description :

Concrete Cube Mould

Technical Specification :

Cube Mould (Cast Iron – for Mortar) 50mm, Single Gang Cube Mould (Cast Iron – for Mortar) 50mm, Three Gang Cube Mould (Mild Steel – for Cement) 70.6mm, with Loose Base Plate Cube Mould (Cast Iron – for Concrete) 75mm Single Gang Cube Mould (Cast Iron – for Concrete) 100mm Light Weight Cube Mould (Cast Iron – for Concrete) 100mm Clamp Type (8.5 KG) Four-Part Cube Mould (Cast Iron – for Concrete) 150mm Light Weight 8KG Cube Mould (Cast Iron – for Concrete) 150mm ISI Marked Cube Mould (Cast Iron -for Concrete) 150mm Clamp Type (18 KG) Four-Part

Cube Mould (Cast Iron – for Concrete) 150mm Clamp Type (16 KG) Two-Part

Apart from Cast Iron and Mild Steel Cube Moulds, we also manufacture Plastic & Polyurethane Cube Moulds.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

