

Product Name :
Concrete Lab Needle Vibrator

Product Code :
NLAB-EDUSCHOOLEQ204012



Description :

Concrete Lab Needle Vibrator

Technical Specification :

The vibrator is mounted on a Spring loaded chair based acting as drip tray and 360° turn table with craving handle for easy portability

Auto air cooled needle

End tip is shouldered with the inner core and crimped for durability

End tip and Inner core made with special steel

Oil cooled rotor provides an even vibration

Specifications:

Size: 25mm X 400 mm, 40 mm X 450 mm and 60 mm X 584 mm

Capacity: 12000 to 15000 cm3

Speed: 2850 / 3000 R.P.M.

Voltage: 440 Volts 3 phase

Power: 3 H.P.

Greaves 4.4 H.P. Air Cooled Diesel Oil Engine complete with Shaft 4 mtrs X 40 mm Needle

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003