

Product Name :
Concrete Lab Permeability Apparatus

Product Code :
NLAB-EDUSCHOOLEQ204005



Description :

Concrete Lab Permeability Apparatus

Technical Specification :

The set-up consist of permeability cell, water reservoir and pressure lines, with glass bottles.

Consists of a metal cylinder with a ledge at the bottom for retaining the specimen, a flange at the top, removable cover plate and a funnel. A rubber gasket is seated in matching grooves, between the cell and the cover plate, to render the joint water tight. Internal working size of the cell will be 100mm dia X 100mm hight.

Water reservoir consists of a length of metal pipe, 50 to 100 mm in diameter and about 500 mm long. The reservoir fitted with a graduated side arm gauge-glass for admitting water and compressed air and a cock for drainage, bleeding and connection to the permeability cell.

Pressure lines heavy duty armoured rubber hose with end fittings.
Complete with glass bottles but without compressor. For specimen size 100 x 100 mm .

Concrete permeability tester for three specimens of 100mm dia x 100mm height or 150mm dia x 150mm height or 150mm cube also available.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003