

Product Name :
Concrete Penetration Test Apparatus

Product Code :
NLAB-TESTINGLABM14002



Description :

Concrete Penetration Test Apparatus

Technical Specification :

Features of Concrete Penetration Test Apparatus

The Apparatus Is Used To Determine The Setting Time Of Concrete By Means Of Penetration

Resistance Measurements On Mortar

Digital Display

Peak Value Automatically Saved

Specifications

Penetration Needle : 100 mm², 50 mm², 20 mm²

Max. Penetration Force : 1000 N

Resolution : 1 N

Power Supply : 220 V AC, 50/60 Hz, Single Phase

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003