

Product Name :
Concrete Railway Sleepers and Bearers Static Testing Machine

Product Code :
NLAB -CONCRLAB120055



Description :

Concrete Railway Sleepers and Bearers Static Testing Machine

Technical Specification :

The 600 kN capacity Flexural Frame is designed to accept all accessories required for flexural, transverse or compression tests.

The load frame provides the stability needed for accurate and repeatable test results over the years of operation.

600 kN capacity Flexural Frame can be connected to any compression machine as a second frame or can be used with any power pack as an independent Flexural Machine.

Flextural test assemblies should be ordered seperately.

Characteristics:

High stability four column design

350 mm piston stroke with safety limit switch

High accuracy load measurement with strain gauge load cells

Can accept wide range of accessories for mentioned standards

The distance between lower rollers can be set up to 1400 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003