

Product Name :
Condenser Aepinus

Product Code :
NLAB-EDUSCHOOLEQ174019



Description :

Condenser Aepinus - Condenser Aepinus, Educational School Laboratory Equipments, Physics Mechanism Instruments, Physics Lab , Physics Mechanics Laboratory Equipment

Technical Specification :

Comprising two metal discs, with variable separation, supported on insulating pillars & provided with insulating handles on the ends of the adjustment slide rods, 4 mm sockets too are drilled onto the ends of the rods.

A rectangular glass plate is supplied which mounts in clips between the two discs & may be easily replaced by a plate of some other material to show the effect of different dielectrics. The whole is mounted on a polished wooden base.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

