

Product Name :
Cone Drop Test Apparatus

Product Code :
NLAB-EDUSCHOOLEQ202002



Description :

Cone Drop Test Apparatus

Technical Specification :

The diameter of the resultant hole is measured with a graduated conical device.

The smaller the hole, the greater the resistance of geosynthetic/ geotextile to damage.

In this test, a geosynthetic or geotextile is clamped to yield a clear diameter of 150 mm and a brass cone of 45 degree angle included, having a weight of 1 kg, is dropped through a height of 500 mm.

When used in combination with other direct tensile test results, it provides a convenient means of qualitative comparison.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003