

**Product Name :**  
Connecting Lanyard Easy Absorb NG 203

**Product Code :**  
NLAB-SAFETYLABP200109



**Description :**

Connecting Lanyard Easy Absorb NG 203

**Technical Specification :**

Length: 1.8 mm.

Strength: 23 kN.

Rope dia: 44 mm.

Twin webbing lanyard with Shock absorber.

With karabiner on one side, and two hooks on the other two sides.

With one side Karabiner and other side two Scaffolding Hooks.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

