

Product Name :
Connecting Lanyard Fall Arrest Single Braided Lanyard NU5001

Product Code :
NLAB-SAFETYLABP20098



Description :

Connecting Lanyard Fall Arrest Single Braided Lanyard NU5001

Technical Specification :

Shock absorber with karabiner and other side scaffolding hooks.

Length: 2 metre from hook to hook.

Rope dia: 10.5 mm.

Strength: 23 kN.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003