

**Product Name :**  
Connecting Lanyard Fall Arrest Single Elasticated Lanyard NU6001

**Product Code :**  
NLAB-SAFETYLABP200102



**Description :**

Connecting Lanyard Fall Arrest Single Elasticated Lanyard NU6001

**Technical Specification :**

Length: Extendable from 1.2 metre to 2 metre.

Shock absorber with karabiner and other side scaffolding hooks.

Extremely lightweight and Compact Lanyard.

Elasticated energy absorbing single Lanyard.

35 mm wide 100% polyester webbing.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003