

Product Name :

Connecting Lanyard Fall Arrest Single Elasticated Lanyard NU6002

Product Code :

NLAB-SAFETYLABP200103

**Description :**

Connecting Lanyard Fall Arrest Single Elasticated Lanyard NU6002

Technical Specification :

Shock absorber with karabiner and other side scaffolding hooks.

Extremely lightweight and Compact Lanyard.

Length: Extendable from 1.2 metre to 2 metre.

Elasticated energy absorbing single Lanyard.

35 mm wide 100% polyester webbing.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003